

CLEANING INSTRUCTIONS FOR PRISMS

DO NOT REMOVE PRISM FROM GLASSES FOR CLEANING

1. Rinse lenses under gentle stream of warm running water.
2. If further cleaning is needed use soft brush (for example a toothbrush) to gently scrub the prism. A little bit of soap can be also be used ensuring adequate rinsing afterwards.
3. Pat or blot dry with a soft, lint-free cloth.

IF PRISM FALLS OFF OR IS REMOVED TO RE-ATTACH

1. Wet spectacle lens and prism with water.
2. Place prism onto lens (the smooth side of the prism should be stuck to the lens).
3. Remove any air bubbles by pressing on the prism gently using a soft tissue or cloth.