

Push-button programme shift

Your hearing instrument has an easy to locate push-button for switching between different programmes.

Programme shift:


Adjust the programme with a quick press.

Your hearing instrument will beep to let you know the programme has changed.

Your hearing instrument can have up to 4 different programmes.


Prog.	No. of beeps	*Prog. assigned by audiologist
1	One	
2	Two	
3	Three	
4	Four	

* To be filled out by your audiologist



To move to the next programme:

Quick press UP



To move back to a previous programme:

Quick Press DOWN